



Point-of-care | Satellite Blood Culture System



Intelligent | Data Management Software



Note: "*" indicates a feature that will be available in the near future.

FA-N Automated Microbial Culture System

Specifications of Culture Bottles

Types and Models
Aerobic culture bottle (AC)
Anaerobic culture bottle (NC)
Pediatric aerobic culture bottle (PC)

Packing Specifications
25 bottles/box, 50 bottles/box

Supported Sample Types
Blood and other sterile body fluids

Blood Sampling Tips

Adults: Each bottle should contain a sample volume of 5–10 mL.
Children: Each bottle should contain a sample volume of 1–3 mL.

Specifications of the Culture System

Normal test environmental conditions
Ambient temperature: 10 °C ~ 30 °C
Relative humidity: 10% ~ 90% RH
Atmospheric pressure: 80.0 kPa ~ 106.0 kPa
Operating altitude: ~400 m ~ 2000 m

References

- [1] Hardy, Liselotte, et al. "Affordable blood culture systems from China: in vitro evaluation for use in resource-limited settings." EBioMedicine 101 (2024).
[2] Aydemir, Özlem, et al. "Comparison of time-to-detection of Mindray TDR and BacT/ALERT® 3D blood culture systems using simulated blood cultures." Acta Clinica Belgica 79.3 (2024): 168-173.



Mindray Medical Official WeChat
Service Center Official WeChat

Service Hotline 400 700 5652
www.mindray.com

P/N: ENG-FA-N Automated Microbial Culture System-210285X6P-20240910

This product brochure and its contents are intended solely for internal communication and reference within the company. For detailed information about the product, please refer to the actual product and relevant documentation. Without written permission from Mindray, no organization or individual may reproduce, quote, distribute, or republish this content in any form. Product specifications are subject to change without prior notice. The latest technical data and test results shall prevail.
©2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

Power Requirements

Power voltage: ~220 V ± 10%
Frequency: 50 Hz ± 1 Hz
Power: 300 VA for FA-N60
600 VA for FA-N120
1200 VA for FA-N240

Dimensions and Weight

Unit	Dimensions (Width x Depth x Height mm)	Weight (kg)
N60 main unit	≤ 540 x 380 x 430	≤ 35
60-bottle extension unit	≤ 540 x 380 x 430	≤ 33
N120/N240 main unit	≤ 560 x 690 x 900	≤ 120
120-bottle extension unit	≤ 290 x 690 x 900	≤ 85



mindray 迈瑞

FA-N Automated Microbial Culture System

Rapid Diagnosis of Bloodstream Infections With Optimized Workflow



mindray 迈瑞
生命科技如此亲近




Enhanced Antibiotic Neutralization and Enriched Culture Media Enable Early Detection of Bloodstream Infections

New Mindray culture bottles for more precise cultivation results

Enhanced Neutralization | For Multiple Antimicrobials

MPC (Multi-Porosity Cross Linked) resin neutralizes 17 categories of antimicrobials.

Examples include *Imipenem*, *Meropenem*, *Piperacillin-tazobactam*, *Cefoperazone-sulbactam*, *Tigecycline*, *Vancomycin*, *Levofloxacin*, *Amikacin*, *Colistin*, *Caspofungin*, *Voriconazole*, and *Amphotericin B*.





Enriched Culture Media | Adequate Microbial Growth Factors

A set (aerobic and anaerobic bottles) is able to recover and detect a broad range of microorganisms (bacteria and fungi)^{[1][2]} in blood and other sterile body fluid samples.

Shorter Time-to-Positivity | Shorter Time-to-Positivity for Slow-Growing Bacteria

Reduced interference from blood cell respiration, along with a customized algorithm for slow-growing bacteria, such as *Talaromyces marneffi*, and *Cryptococcus neoformans*.

A technological breakthrough: development of an algorithm for slow-growing bacteria



Shorter time-to-positivity



An Integrated and Intelligent POC Information System Enables Effective and Qualified Results

Seamless Connection for Faster Result Delivery

Integrated | In-time Blood Volume Monitoring Technology

